

University of Groningen

Exploiting Catalytic Promiscuity for Biocatalysis

Miao, Yufeng

IMPORTANT NOTE: You are advised to consult the publisher's version (publisher's PDF) if you wish to cite from it. Please check the document version below.

Document Version

Publisher's PDF, also known as Version of record

Publication date:

2015

[Link to publication in University of Groningen/UMCG research database](#)

Citation for published version (APA):

Miao, Y. (2015). *Exploiting Catalytic Promiscuity for Biocatalysis: Carbon-Carbon Bond Formation by a Proline-Based Tautomerase*. [Thesis fully internal (DIV), University of Groningen]. University of Groningen.

Copyright

Other than for strictly personal use, it is not permitted to download or to forward/distribute the text or part of it without the consent of the author(s) and/or copyright holder(s), unless the work is under an open content license (like Creative Commons).

The publication may also be distributed here under the terms of Article 25fa of the Dutch Copyright Act, indicated by the "Taverne" license. More information can be found on the University of Groningen website: <https://www.rug.nl/library/open-access/self-archiving-pure/taverne-amendment>.

Take-down policy

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

Downloaded from the University of Groningen/UMCG research database (Pure): <http://www.rug.nl/research/portal>. For technical reasons the number of authors shown on this cover page is limited to 10 maximum.

4

Practical method: Asymmetric Michael-type Additions of Acetaldehyde to Nitroolefins Catalyzed by 4-Oxalocrotonate Tautomerase (4-OT) Yielding Valuable γ -Nitroaldehydes


Edzard M. Geertsema, Yufeng Miao, and Gerrit J. Poelarends

In Practical Methods in Biocatalysis and Biotransformations – volume 3 (Eds.: J. Whittal, P. Sutton, and W. Kroutil), 2015, Wiley, UK, in press.



Introduction

γ -Nitroaldehydes are useful precursors for γ -aminobutyric acid (GABA) derivatives, several of which represent marketed pharmaceuticals.¹ The Michael-type addition of acetaldehyde to nitroolefins gives convenient access towards relevant chiral γ -nitroaldehydes (Scheme 1). The first successful methodology for this type of reaction was developed in the field of organocatalysis^{2,3} and documented in the literature in 2008. We lately presented a biocatalytic procedure for the Michael-type addition of acetaldehyde to various nitroolefins, which involves the homohexameric enzyme 4-oxalocrotonate tautomerase (4-OT).^{4,5,6} 4-OT is a member of the tautomerase superfamily and is part of a catabolic pathway for aromatic hydrocarbons in the soil bacterium *Pseudomonas putida* mt-2. 4-OT is a small enzyme (only 62 amino acid residues per monomer) and is characterized by a catalytic N-terminal proline residue (Pro-1) that resides in the active site.^{7,8,9,10} We here describe the production and purification of the enzyme 4-OT, as well as its use as a catalyst in the Michael-type addition of acetaldehyde to nitroolefins **1a-b** to produce γ -nitroaldehydes **2a-b** (Scheme 1).^{4,5,6}



Scheme 1. Michael-type addition of acetaldehyde to nitroolefins **1a-b** yielding chiral γ -nitroaldehydes **2a-b**.
* = chiral center.

Materials and equipment for expression and purification of 4-OT

- (Conical-bottom) centrifuge tube (50 mL, CELLSTAR)
- Ammonium sulphate (3.2 M), 50 mL
- Ampicillin (100 mg/ml stock solution in water, filter sterilized)
- Centrifuge capable of reaching 15000 rpm with temperature control
- DEAE-sepharose (CL-6B, Sigma-Aldrich) column (10×1.0 cm filled with ~8 ml DEAE-sepharose resin)
- Dialysis membrane with MWCO 2000 Da (Spectrum labs)
- *E. coli* BL21(DE3) containing pET20b(4-OT) expression vector (frozen glycerol stock at -80°C)

- Erlenmeyer flask (3 L)
- Eppendorf tubes (1.5/2 mL)
- LB medium (1 L, sterilized)
- Magnetic stirrer
- PD-10 sephadex G-25 gel filtration column (GE Healthcare Bio-sciences AB)
- Phenyl-sepharose (CL-4B, Sigma-Aldrich) column (10×1.0 cm filled with ~8 ml phenyl-sepharose resin)
- Quartz cuvette with 1 mm or 10 mm light path (Hellma Analytics)
- Rotary shaker with temperature control
- Sodium dihydrogenphosphate (NaH_2PO_4 , 20 mM, pH 7.3), 100 mL
- Sodium dihydrogenphosphate (NaH_2PO_4 , 10 mM, pH 8.0) containing ammonium sulphate [$(\text{NH}_4)_2\text{SO}_4$, 1.6 M], 100 mL
- Sodium dihydrogenphosphate [NaH_2PO_4 , 10 mM, pH 8.0] containing sodium sulphate [Na_2SO_4 , 0.5 M], 100 mL
- Sodium dihydrogenphosphate [NaH_2PO_4 , 10 mM, pH 8.0], (buffer A), 250 mL
- Sodium dodecyl sulfate (SDS) gels containing polyacrylamide (10%)
- Sterile loop
- Syringe filter, pore size 0.2 μm (FP30/0.2, Whatman)
- UV-VIS spectrophotometer
- Vortex mixer

Procedure for expression and purification of 4-OT

1. Using a sterile loop, cells from a frozen glycerol stock of *E. coli* BL21(DE3) harboring the pET20b(4-OT) expression vector were used to inoculate 5 mL LB medium, in a 50 mL conical-bottom centrifuge tube, containing ampicillin (100 $\mu\text{g}/\text{mL}$).
2. After overnight growth with shaking (250 rpm) at 37°C, this culture was used to inoculate the rest of LB medium (995 mL) containing ampicillin (100 $\mu\text{g}/\text{mL}$) in a 3 L Erlenmeyer flask. The culture was incubated overnight at 37°C with vigorous shaking (200 rpm) on a rotary shaker.
3. The culture was divided into two equal portions and cells were harvested by centrifugation (20 min at 4500 rpm) after an OD_{600} value of ~4.5 was reached.¹¹
4. Cells harvested from 0.5 L of culture were resuspended in ~10 mL of buffer A by using a vortex mixer. [Optional: for storing the cells at -20°C, leave the cells in 10 mM NaH_2PO_4 buffer (pH 8.0, buffer A).]
5. The cell suspension was transferred to sonication tubes (centrifuge tubes) and sonicated at high power output for 8-10 min to disrupt the cells. After centrifugation (45 min at 15000 rpm), the soluble fraction (supernatant) containing 4-OT was separated from the insoluble part (pellet) [Note: Keep the cell suspension in a cooling bath (e.g. water/ice) during the sonication process to avoid heat inactivation of 4-OT.]

6. The supernatant that contains 4-OT was then loaded onto a DEAE-sepharose column, which had been previously equilibrated with buffer A, using gravity flow. The flow-through was discarded.
7. The column was first washed with 3×10 mL buffer A and then the protein was eluted by gravity flow using buffer A containing 0.5 M Na₂SO₄ (12 mL). Fractions (1.5 mL × 8) from the elution step were collected in Eppendorf tubes (1.5 mL) and the presence of 4-OT in each fraction was examined by SDS-PAGE.
8. The fractions containing 4-OT were then combined and mixed with (NH₄)₂SO₄ (final concentration: 1.6 M). After gentle stirring with a magnetic stir bar at 4°C for 2 h, the precipitate was pelleted by centrifugation (20 min at 13300 rpm, 4°C).
9. The supernatant was filtered (syringe filter with 0.2 μm pore size) and loaded onto a phenyl-sepharose column that had been previously equilibrated with buffer A containing 1.6 M (NH₄)₂SO₄.
10. The supernatant was left to flow through the column and the flow-through was discarded. The column was then washed with 3×10 mL of loading buffer [buffer A containing 1.6 M (NH₄)₂SO₄], after which the protein was eluted by gravity flow using buffer A (12 mL). Fractions of the wash step (~30 mL) and the elution step (~1.5 mL × 8) were analyzed by SDS-PAGE. [Note: wild-type 4-OT has little interaction with the phenyl-sepharose column and should therefore elute as homogenous protein (>95% purity as assessed by SDS-PAGE) in the wash step.]
11. Fractions that contained pure 4-OT were combined and the buffer was exchanged against 20 mM NaH₂PO₄ buffer (pH 7.3) using a pre-packed PD-10 sephadex G-25 gel filtration column. [Optional: An additional dialysis step is recommended to reduce the high salt concentration (1.6 M (NH₄)₂SO₄) in the sample before exchanging the buffer on a PD-10 column. In a typical dialysis procedure, the fractions containing pure 4-OT were combined and transferred into a dialysis membrane (MWCO 2000 Da). The dialysis membrane containing the sample was then placed in a 5 L beaker containing ~4 L 20 mM NaH₂PO₄ buffer (pH 7.3) with gentle stirring on a magnetic stirrer overnight.]
12. The concentration of purified 4-OT in 20 mM NaH₂PO₄ buffer (pH 7.3) was determined with a spectrophotometer using the method of Waddell.¹² A typical yield of purified 4-OT that can be isolated from 0.5 L of cell culture following this protocol is 25–50 mg.
13. The purified protein was divided into aliquots (100–200 μl) and stored at 4°C until further use. [Optional: Purified 4-OT can be stored for at least one month upon freezing the protein with liquid nitrogen and when kept at -80°C.]

Materials and equipment for synthesis of 2a

- *Trans*-β-nitrostyrene **1a** (18.0 mg, commercially available, Sigma-Aldrich)
- Acetaldehyde (132 mg, ~0.17 mL)
- EtOH (6 mL)
- 20 mM NaH₂PO₄ buffer, pH 5.5 (~54.0 mL)

- 4-OT (11.3 mg from a stock solution in 20 mM NaH_2PO_4 buffer, pH 7.3)
- Diethyl ether (120 mL)
- MgSO_4 (~0.5 g)
- Deuterated chloroform (CDCl_3 , 0.65 mL for ^1H NMR analysis)
- Balance
- Measuring cylinder (50 mL)
- Glass flask (150 mL)
- Standard pipette for 4-OT addition
- 1 mL Syringe (NORM-JECT) for acetaldehyde addition
- UV-VIS spectrophotometer
- Cuvette (preferably path length = 1 mm, volume = 300 μL . See paragraph 'analytical methods')
- Two vivaspin columns (Sartorius Stedim Biotech S.A., France) with a 5000 Da molecular weight cut-off filter
- Separatory funnel (250 mL)
- Two Erlenmeyer flasks (100 mL and 250 mL)
- Magnetic stirrer
- Magnetic stir bar (for drying with MgSO_4 during work up procedure)
- Funnel with filter paper
- Round bottom glass flask (250 mL)
- Rotary evaporator connected to a vacuum pump
- NMR tube
- NMR spectrometer

Procedure for synthesis of **2a** (from reference 5)

1. In a glass flask (150 mL), a solution of nitroolefin **1a** (18.0 mg, 0.12 mmol) in EtOH (6.0 mL) was added to a mixture of acetaldehyde (132 mg, 3.0 mmol) and 4-OT (11.3 mg from a stock solution in 20 mM NaH_2PO_4 buffer pH 7.3, 1.7×10^{-3} mmol, 1.4 mol%) in 20 mM NaH_2PO_4 buffer, pH 5.5 (final volume of reaction mixture: 60 mL. Final ratio water/EtOH: 90/10 v/v).
2. The mixture was incubated at room temperature and reaction progress was monitored by recording UV-spectra of aliquots taken from the reaction mixture after regular time intervals. Note: to avoid inactivation of the enzyme, it is recommended not to stir the reaction mixture with a stirring magnet but rather to gently shake the mixture every 30 min.
3. After 2.0 h, all **1a** was converted. The reaction mixture was divided into two equal portions, transferred to two Vivaspin columns (Sartorius Stedim Biotech S.A., France) with a 5000 Da molecular weight cut-off filter, and centrifuged (4000 rpm, RT) to remove 4-OT.
4. The flow-throughs were collected, combined, and extracted with diethyl ether (3 \times 40 mL). The combined organic layers were dried with MgSO_4 , filtered, and concentrated

in vacuo with a rotary evaporator to yield **2a** with high purity (16.2 mg, 8.4×10^{-2} mmol, 70%) as a colorless oil.

5. The ^1H NMR spectroscopic data of **2a** were in agreement with published data.^{13,14} Enantiomeric excess was determined by derivatization of the aldehyde moiety of **2a** into the corresponding cyclic acetal.⁴ Normal phase HPLC analysis of derivatized **2a** using a Chiracel OD column revealed an enantiomeric excess of 81% in favor of the (*S*)-enantiomer. HPLC parameters: eluent: (*n*-heptane/*i*-PrOH 90:10, 40°C), flow rate: 1 mL/min, UV detection at 210 nm, t_{R} : (minor) = 13.0 min, (major) = 15.4 min. The HPLC data are in accordance with the literature.¹⁴

Materials and equipment for synthesis of **2b**

- Nitroolefin **1b**: (*E*)-2-(cyclopentyloxy)-1-methoxy-4-(2-nitrovinyl)benzene (31.6 mg)¹⁵
- Acetaldehyde (132 mg, ~0.17 mL)
- DMSO (24 mL)
- 20 mM NaH_2PO_4 buffer, pH 5.5 (~36.0 mL)
- 4-OT (30 mg from stock solution in 20 mM NaH_2PO_4 buffer, pH 7.3)
- Water (18 mL)
- Chloroform (48 mL)
- MgSO_4 (~0.5 g)
- Deuterated chloroform (CDCl_3 , 0.65 mL for ^1H NMR analysis)
- Balance
- Measuring cylinder (50 mL)
- Glass flask (150 mL)
- Standard pipette for 4-OT addition
- 1 mL Syringe (NORM-JECT) for acetaldehyde addition
- UV-VIS spectrophotometer
- Cuvette (preferably path length = 1 mm, volume = 300 μL . See paragraph 'analytical methods')
- Six CELLSTAR tubes (15 mL, polypropylene (PP), conical bottom)
- Acid-resistant CentriVap vacuum concentrator (Labconco, 78100 series) connected to a cold trap and a vacuum pump
- Separatory funnel (100 mL)
- Two Erlenmeyer flasks (25 mL and 100 mL)
- Magnetic stirrer
- Magnetic stir bar (for drying with MgSO_4 during work up procedure)
- Funnel with filter paper
- Round bottom glass flask (100 mL)
- Rotary evaporator connected to a vacuum pump
- NMR tube
- NMR spectrometer

Procedure for synthesis of **2b** (from reference 6)

1. In a glass flask (150 mL), a solution of nitroolefin **1b** (31.6 mg, 0.12 mmol) in DMSO (6.0 mL) was added to a mixture of acetaldehyde (132 mg, 3.0 mmol) and 4-OT (30 mg from a stock solution in 20 mM NaH₂PO₄ buffer pH 7.3, 4.4×10⁻³ mmol, 3.7 mol%) in 20 mM NaH₂PO₄ buffer (pH 5.5) and DMSO (final volume of reaction mixture: 60.0 mL. Final ratio water/DMSO: 60/40 v/v).
2. Note: DMSO was required as a co-solvent to achieve sufficient solubility of nitroolefin **1b** and product **2b** in aqueous buffer. DMSO enhances solubility of **1b** and **2b**, does not impede catalytic activity of 4-OT, and does not chemically react with any of the reagents (see reference 6 for full explanation).
3. The mixture was incubated at room temperature and reaction progress was monitored by recording UV-spectra of aliquots taken from the reaction mixture after regular time intervals. Note: to avoid inactivation of the enzyme, it is recommended not to stir the reaction mixture with a stirring magnet but rather to gently shake the mixture every 30 min.
4. After 2.5 h, all **1b** was converted and the reaction mixture was divided over six polypropylene tubes (15 mL, CELLSTAR). The solvents were evaporated using an acid-resistant CentriVap vacuum concentrator (55°C, overnight).
5. Water (3 mL) and chloroform (3 mL) were added to each tube. The dry residues were dissolved/suspended by vigorous stirring (additional scraping with a spatula may be required). The combined water and chloroform layers were separated in a separatory funnel. The water layer was extracted with chloroform (3×10 mL). The combined organic layers were dried with MgSO₄, filtered, and concentrated *in vacuo* with a rotary evaporator to yield **2b** with high purity (23.6 mg, 7.7×10⁻² mmol, 64%) as a colorless oil.
6. The ¹H NMR spectroscopic data of **2b** were in agreement with published data.² Enantiomeric excess was determined by derivatization of the aldehyde moiety of **2b** into its corresponding methyl ester.¹⁶ Normal phase HPLC analysis of derivatized **2b** using a Chiralpak IB column revealed an enantiomeric excess of 96% in favor of the (*S*)-enantiomer.¹⁶ HPLC parameters: eluent: (*n*-heptane/*i*-PrOH 95:5, 25°C), flow rate: 1 mL/min, UV detection at 220 nm, t_R: (minor) = 32.2 min, (major) = 34.9 min.

Analytical methods

The biocatalytic conversions of nitroolefins **1a-b** into **2a-b** were monitored by UV spectroscopy. After regular time intervals (typically 30 min), 300 μ L of reaction mixture was transferred into a 300 μ L cuvette with a path length of 1 mm (see note below) and a UV-spectrum was recorded (200–500 nm). A decrease of absorbance at λ_{\max} of **1a-b** indicates the depletion of **1a-b** ($\lambda_{\max,1a} = 320$ nm, $\epsilon_{\max} = 14.4$ mM⁻¹ cm⁻¹; $\lambda_{\max,1b} = 378$ nm, $\epsilon_{\max} = 13.5$ mM⁻¹ cm⁻¹). After recording the UV-spectrum, the contents of the cuvette were recombined with the reaction mixture (with a syringe). The work up procedure was initiated when the absorbance at λ_{\max} of the nitroolefin (**1a-b**) had totally vanished which indicates full conversion of the nitroolefin (**1a-b**).

Note: obviously one may use a cuvette with a path length of 10 mm and volume of ~1 mL instead. However, one should then take an aliquot of 100 μ L of the reaction mixture and dilute with 900 μ L of buffer (20 mM NaH_2PO_4 , pH 5.5).

Results and conclusion

The enzyme 4-OT can be readily produced and purified. About 50 to 100 mg of 4-OT can be obtained from 1 L of cell culture. The enzyme can be stored at -80°C in concentrations of 10 to 20 mg/mL (in 20 mM NaH_2PO_4 buffer, pH 7.3) for at least one month without significant loss of catalytic activity. 4-OT can be easily tested on whether it is fully catalytically active by assaying its natural tautomerase activity with the substrates 2-hydroxyomuconate¹⁷ or phenyl(enol)pyruvate¹⁸ and comparing with literature data.

The enzyme 4-OT promiscuously catalyzes the Michael-type addition of acetaldehyde to nitroolefins **1a-b** yielding γ -nitroaldehydes **2a-b** (Scheme 1). The methodology is characterized by relatively low catalyst loading (compared to 10–20 mol% in organocatalysis),^{2,3,19,20,21,22} high stereoselectivities, and aqueous reaction media (Table 1). Reaction times of 2.5 h or less were sufficient for full conversion of nitroolefins **1a-b**, respectively. Side products were hardly formed as revealed by straightforward work-up procedures without column chromatography which yielded **2a-b** with high purity (¹H NMR spectroscopy). Respective yields of 70% and 64% were established for **2a-b** while e.e. values of 81% and 96%, respectively, were found. Catalysis of the Michael-type additions takes place in the active site of 4-OT and its characteristic N-terminal proline is a key catalytic residue as demonstrated by various control experiments.^{4,5,6} 4-OT fully retains its ‘Michaelase’ activity in aqueous solvent systems containing up to 50% DMSO (v/v).⁶ This methodology is therefore not only restricted to water-soluble substrates but also allows application of poorly water-soluble chemicals such as **1b** since its solubility in aqueous solvent systems is greatly enhanced in the presence of DMSO. In addition to acetaldehyde and aromatic nitroolefins, 4-OT also accepts linear aldehydes up to octanal as donors⁵ and aliphatic nitroolefins as acceptors⁶ for Michael-type addition reactions.

Table 1. Conditions for, and results of, 4-OT-catalyzed Michael-type additions of acetaldehyde to nitroolefins **1a** and **1b** in 20 mM NaH_2PO_4 buffer (pH 5.5) and co-solvent yielding **2a** and **2b** (Scheme 1).

| nitroolefin | 4-OT (mol%) | co-solvent (v/v) | reaction time (h) | product | yield (%) | e.e. (%) |
|-------------|-------------|------------------|-------------------|-----------|-----------|----------|
| 1a | 1.4 | EtOH 10% | 2 | 2a | 70 | 81 |
| 1b | 3.7 | DMSO 40% | 2.5 | 2b | 64 | 96 |

References

1. Maltsev, O. V.; Kucherenko, A. S.; Beletskaya, I. P.; Tartakovsky, V. A. and Zlotin, S. G. *Eur. J. Org. Chem.* **2010**, 2927-2933 and references cited herein.
2. García-García, P.; Ladépêche, A.; Halder, R. and List, B. *Angew. Chem. Int. Ed.* **2008**, *47*, 4719-4721. Reagents: acetaldehyde and **1a**; Mol% catalyst: 20%; reaction time: 15 h; yield: 51%; e.e.: 92%
3. Hayashi, Y.; Itoh, T.; Ohkubo, M. and Ishikawa, H. *Angew. Chem. Int. Ed.* **2008**, *47*, 4722-4724. Reagents: acetaldehyde and **1a**; Mol% catalyst: 10%; reaction time: 18 h; yield: 75%; e.e.: 96%
4. Zandvoort, E.; Geertsema, E. M.; Baas, B. J.; Quax, W. J. and Poelarends, G. J. *Angew. Chem. Int. Ed.* **2012**, *51*, 1240-1243.
5. Miao, Y.; Geertsema, E. M.; Tepper, P. G.; Zandvoort, E. and Poelarends, G. J. *ChemBioChem* **2013**, *14*, 191-194.
6. Geertsema, E. M.; Miao, Y.; Tepper, P. G.; De Haan, P.; Zandvoort, E. and Poelarends, G. J. *Chem. Eur. J.* **2013**, *19*, 14407-14410.
7. Whitman, C. P. *Arch. Biochem. Biophys.* **2002**, *402*, 1-13.
8. Poelarends, G. J.; Puthan Veetil, V. and Whitman, C. P. *Cell. Mol. Life Sci.* **2008**, *65*, 3606-3618.
9. Poelarends, G. J. and Whitman, C. P. *Bioorg. Chem.* **2004**, *32*, 376-392.
10. Baas, B. J.; Zandvoort, E.; Geertsema, E. M. and Poelarends, G. J. *ChemBioChem* **2013**, *14*, 917-926.
11. Sambrook, J.; Fritsch, E. F. and Maniatis T. in *Molecular Cloning: A Laboratory Manual* (2nd ed.), Cold Spring Harbor Laboratory Press, **1989**.
12. Waddell, W. J. *J. Lab. Clin. Med.* **1956**, *48*, 311-314.
13. Gotoh, H.; Ishikawa, H. and Hayashi, Y. *Org. Lett.* **2007**, *9*, 5307-5309.
14. Mager I. and Zeitler K. *Org. Lett.* **2010**, *12*, 1480-1483.
15. Hynes, P. S.; Stuppel, P. A. and Dixon, D. J. *Org. Lett.* **2008**, *10*, 1389-1391.
16. Palomo, C.; Landa, A.; Mielgo, A.; Oiarbide, M.; Puente, A. and Vera, S. *Angew. Chem. Int. Ed.* **2007**, *46*, 8431-8435.
17. Whitman, C. P.; Aird, B. A.; Gillespie, W. R. and Stolowich, N. J. *J. Am. Chem. Soc.* **1991**, *113*, 3154-3162.
18. Johnson, W. H.; Wang, S. C.; Stanley, T. M.; Czerwinski, R. M.; Almrud, J. J.; Poelarends, G. J.; Murzin, A. G. and Whitman, C. P. *Biochemistry* **2004**, *43*, 10490-10501.
19. Alza, E. and Pericas, M. A. *Adv. Synth. Catal.* **2009**, *351*, 3051-3056; Reagents: acetaldehyde and **1a**; Mol% catalyst: 10%; reaction time: 72 h; yield: 55%; e.e.: 96%.
20. Riente, P.; Mendoza, C. and Pericas, M. A. *J. Mater. Chem.* **2011**, *21*, 7350-7355; Reagents: acetaldehyde and **1a**; Mol% catalyst: 10%; reaction time: 48 h; yield: 58%; e.e.: 0%.
21. Jentzsch, K. I.; Min, T.; Etcheson, J. I.; Fettingter, J. C. and Franz, A. K. *J. Org. Chem.* **2011**, *76*, 7065-7075; Reagents: acetaldehyde and **1a**; Mol% catalyst: 20%; reaction time: 3 h; yield: 75%; e.e.: 95%.
22. Qiao, Y.; He, J.; Ni, B. and Headley, A. D. *Adv. Synth. Catal.* **2012**, *354*, 2849-2853. Reagents: acetaldehyde and **1a**; Mol% catalyst: 10%; reaction time: 22 h; yield: 57%; e.e.: 91%.



