

University of Groningen

STAT5 is required for long-term maintenance of normal and leukemic human stem/progenitor cells

Schepers, Hein; van Gosliga, Djoke; Wierenga, Albertus T. J.; Eggen, Bart J. L.; Schuringa, Jan Jacob; Vellenga, Edo

Published in:
Blood

DOI:
[10.1182/blood-2006-08-039073](https://doi.org/10.1182/blood-2006-08-039073)

IMPORTANT NOTE: You are advised to consult the publisher's version (publisher's PDF) if you wish to cite from it. Please check the document version below.

Document Version
Publisher's PDF, also known as Version of record

Publication date:
2007

[Link to publication in University of Groningen/UMCG research database](#)

Citation for published version (APA):

Schepers, H., van Gosliga, D., Wierenga, A. T. J., Eggen, B. J. L., Schuringa, J. J., & Vellenga, E. (2007). STAT5 is required for long-term maintenance of normal and leukemic human stem/progenitor cells. *Blood*, 110(8), 2880-2888. <https://doi.org/10.1182/blood-2006-08-039073>

Copyright

Other than for strictly personal use, it is not permitted to download or to forward/distribute the text or part of it without the consent of the author(s) and/or copyright holder(s), unless the work is under an open content license (like Creative Commons).

Take-down policy

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

Downloaded from the University of Groningen/UMCG research database (Pure): <http://www.rug.nl/research/portal>. For technical reasons the number of authors shown on this cover page is limited to 10 maximum.

Day 1

- Isolate 1.5×10^6 AML CD34⁺ cells with MACS column or sort by FACS and incubate (0.5×10^6 /ml) in RPMI1640 supplemented with 10% FCS, IL-3, G-CSF and TPO (each 20 ng/ml) for 4 hours, 37°C, 5% CO₂
- Coat 3 wells from a 12-wells-plate with retronectin: Dilute with PBS to 50 µg/ml, 0.5 ml per well of a non-tissue culture 12-well plate. Incubate at room temperature for 2 hrs, then remove retronectin and block with 2% BSA/PBS for 30 min. Wash plate twice with PBS and keep in fridge until use.
- After 4 hours split AML CD34⁺ cells in 3 groups (at least 1.5×10^5 cells per group, but not more than 5×10^5 in 500 µl) and apply each group to one coated well:

Group 1: No virus control (for setting FACS gates)

Group 2: Transduce with Renilla RNAi control vector

Group 3: Transduce with STAT5 RNAi vector

- Add 500 µl HPGM to group 1
Add 500 µl Renilla RNAi lentivirus supernatant to group 2
Add 500 µl STAT5 RNAi lentivirus supernatant to group 3
- Add to each well to final volume/concentration:
10% FCS
20 ng/ml hIL-3
20 ng/ml TPO
20 ng/ml GCSF
4 µg/ml polybrene
- Incubate overnight at 37°C, 5% CO₂ (round 1)

Furthermore:

Culture MS5 stromal cells in T75 in α -MEM supplemented with 10% FCS, so that 80% confluence is reached at day 2

Day 2

- Repeat transduction procedure in morning (round 2). Make virus supernatants (or HPGM) with 10% FCS and growth factors (each 20 ng/ml) and polybrene (4 ng/ml) (again in 500 µl) as above.
- Add 500 µl to each well (Washing of cells is not necessary).
- Incubate for 8 hrs
- Repeat this procedure in the evening (round 3)

Furthermore:

- Coat 3x T25 with 0.1 % gelatin for 2 hrs at room temperature
- Trypsinize T75 with MS5, resuspend in 15 ml α -MEM with 10% FCS and plate 5 ml in gelatin-coated T25s (these MS5 cultures will reach confluency on day 3)

Day 3

- Wash AML CD34⁺ cells from each group 3 times with PBS and resuspend in 1.5 ml of Long Term Culture (LTC) medium (α MEM supplemented with heat-inactivated 12.5% FCS, heat-inactivated 12.5% Horse serum, penicillin and streptomycin, 200 mM Glutamine, 57.2 μ M β -mercaptoethanol and 1 μ M hydrocortisone) supplemented with IL-3, G-CSF and TPO (each 20 ng/ml).
- Take a 50 μ l aliquot for FACS to assess transduction efficiency (use group 1 cells to set the gates).
- Plate equal amounts of AML CD34⁺ cells on MS5 in LTC medium (at least 1.5×10^5 cells per T25, no more than 5×10^5). Sorting of YFP/GFP⁺ cells is not necessary. The untransduced GFP/YFP⁻ cells within each culture serve as an internal control for each culture.

Every week

- Remove weekly half of the suspension (2.5 ml) from the cultures and add 2.5 ml fresh LTC medium supplemented with IL-3, G-CSF and TPO (each 20 ng/ml) (demipopulate). The removed half can be used for FACS analysis, cell count, cytopins etc.

Addendum:

1. The GFP/YFP percentage was usually determined after 2 days expansion on MS-5 stromal layer rather than directly after transduction. Then, GFP/YFP expression will have increased until steady levels and a small sample can be analyzed for proper transduction efficiency.
2. Depending upon the expansion and growth of the AML cells, the MS5 stromal cell layer may be exhausted after 4 to 5 weeks of culture. By isolating all the cells from the culture (including the trypsinized adherent fraction, consisting of MS5 and AML CAFs) and replating 10% of these cells into a new T25 (freshly coated with new MS5 cells), new cultures can be initiated and make the investigation of long term AML cultures possible.